



Mercury 2V dry heat sterilizer

- 1. Modern design and ease of use:** the dry heat sterilizer opens from the top like a drawer, allowing easy access to the instruments. All controls are intuitive, with clear indicators showing temperature and time.
- 2. Timer with alert function:** when sterilization is complete, in addition to the visual display on the indicator, the dry heat sterilizer will beep to indicate the completion of the process.
- 3. Effective sterilization modes:** Mercury 2V has three automatic sterilization modes (160°C, 180°C, 200°C). Each of them has an appropriate time to guarantee safe cleaning of instruments.
- 4. High safety standards:** the dry heat sterilizer has been tested and certified in Ukrainian and European laboratories, which confirms its high quality and reliability.
- 5. Intelligent safety systems:** automatic shutdown when the temperature in the sterilization chamber or on the outer surface of the lid is exceeded provides an additional level of protection. The emergency opening system allows you to open the lid of the dry heat sterilizer with an electric lock, even without electrical power.

Dimensions of the dry heat sterilizer, mm	274 x 228 x 300
Sterilization chamber dimensions, mm	110x115x245
The type of loading	вертикальный
Temperature deviation from the set value, °C	-1/*5°
Drying modes	160°C, 180°C, 200°C
Heating time to sterilization temperature, min	8-10
Maximum power consumption, kW	0,4
Weight of the dryer, kg	8.6
Capacity of kraft bags	3
Emergency shutdown	Yes
The presence of a fan	Yes
Warranty	24 months

1. The **ÜLKA Mercury 2V** dry heat sterilizer is distinguished by its innovative design, where the **buttonless door opens like a magic box** where sterilized instruments become a real treasure. This feature allows for **easy access to the sterilization chamber**, making it a real pleasure to work with.
2. Thanks to its **compact size**, the Mercury 2V dry heat sterilizer offers great convenience as it takes up very little space. This makes it easy to move around, **making it ideal for small studios or private salons**.
3. The dryer is equipped with a **timer and an operating indicator** that sounds when the process is complete, making the use of the device intuitive and convenient. The device is also equipped with three automatic sterilization modes (**160°C, 180°C, 200°C**), each of which meets sanitary standards and has a set time interval.
4. Also, taking into account the system of opening the hot air sterilizer from the top, we added a stopper for the opening lid - this helps to protect the sterilizer's body from damage and scratches. We have developed a **unique technology** for opening the dry heat sterilizer's cabinet by adding an **electric lock** - now, the nail artist does not need to exert much effort by pressing the heavy dry heat lid - it opens easily by pressing a button on the body. In addition, we have provided an **emergency lid opening mechanism** - even in the absence of light, the dry heat cabinet will open at the touch of a button, and the nail technician can take out the sterilized tools.



Your reliable assistant in sterilizing instruments! The Mercury 2V is an effective solution for professionals that guarantees the safety of the nail technician and the client.